

Catalog Number: IRT002

Anti-His tag antibody


Package size: 200, 100, 25 µl
Store at: -20°C
MW (kDa):

Overview

Product Name	His tag antibody
Product Number	IRT002
Clonality	Polyclonal
Host	Rabbit
Immunogen	KLH-conjugated HHHHHH (6X His)
Recommended Applications Dilutions	Western Blot 1:2000-1:4000 Immunofluorescence 1:500 – 1:1000
Storage Buffer	100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1.01 mg/ml
Purity	Affinity column purified
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

Image

IRT002 anti His-tag antibody WB image



His-tag antibody detects His-tagged protein by western blot analysis. The membrane was then blocked for an hour before being incubated with IRT002 overnight at 4°C. His-tag antibody (IRT002) diluted at 1:2000.

IReal Biotechnology Co., Ltd.