

Catalog Number: IR306-909-PE


Anti L-selectin / CD62L (**PE conjugated**) Antibody

Package size: 25, 100 μ l
Store at: 4°C

Overview

| | |
|------------------------------------|---|
| Product Name | L-selectin / CD62L (PE conjugated) antibody |
| Product Number | IR306-909-PE |
| Gene Description | selectin, lymphocyte |
| Clonality | Polyclonal |
| Host | Rabbit |
| Species Reactivity | Human, mouse, rat |
| Recommended Applications Dilutions | Immunofluorescence 1:200 – 1:300 Flow Cytometry 1:50 – 1:100 |
| Storage Buffer | PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN3 was added as a preservative |
| Concentration | 1 mg/ml |
| Purity | Affinity column purified |
| Storage | This product is stable for 12 months when stored at 4°C. Do not aliquot the antibody. Protect from light. Do not freeze. |
| Notes | Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user. |

L-selectin / CD62L **PE conjugated** antibody FACS image



Product: L-selectin / CD62L **PE conjugated** antibody
Catalog: IR306-909-PE
Cell: PBMC
L-selectin / CD62L antibody: 1-100
Data from: POLOWANG BIOTECH INC

iReal Biotechnology Co., Ltd.