

Catalog Number: IR362

Anti-Procalcitonin / PCT antibody

Package size: 25, 100 μ l
Store at: -20°C
MW (kDa): -

Overview

Product Name	Procalcitonin / PCT antibody
Product Number	IR362
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human
Recommended Applications Dilutions	Western Blot 1:500 – 1:1000 Immunofluorescence 1:200 – 1:400
Storage Buffer	100mM Tris Glycine, 30% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1.14 mg/ml
Purity	Affinity column purified
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IR362 anti-Procalcitonin / PCT antibody WB image

IR362
Procalcitonin / PCT antibody



50 10 2 0.2 (ng)

Anti-Procalcitonin / PCT antibody at 1/1000 dilution

Specificity and sensitivity of **Procalcitonin / PCT** Antibody (IR362)

Nitrocellulose membrane was then blocked for an hour before being incubated with IR362 overnight at 4°C.

iReal Biotechnology Co., Ltd.

