

Catalog Number: IRCP2312M13


Anti M13 Phage antibody

Package size: 25, 100 μ l
Store at: -20°C
MW (kDa): -

Overview

Product Name	M13 Phage antibody
Product Number	IRCP2312M13
Gene Description	Capsid protein G8P (M13)
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Virus
Recommended Applications Dilutions	Western Blot 1:2000 – 1:4000 ELISA 1:5000 – 1:20000
Storage Buffer	100mM Tris Glycine, 30% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRCP2312M13 anti-M13 Phage antibody WB image



Sensitivity analysis: Coating the indicated amounts of M13 bacteriophages and stained with IRCP2312M13 anti-M13 Phage antibody.

IReal Biotechnology Co., Ltd.