

Catalog Number: IRM056

Anti-human CD45 [HI30] monoclonal antibody


Package size: 25, 100 μ l
Store at: -20°C
MW (kDa):

Overview

Product Name	human CD45 [HI30] monoclonal antibody
Product Number	IRM056
Clonality	Monoclonal
Host	Mouse (mouse IgG1)
Species Reactivity	Human
Recommended Applications Dilutions	Immunofluorescence 1:100 – 1:200 Immunohistochemistry (Paraffin) 1:100
Storage Buffer	Phosphate-buffered solution, pH 7.2, 50% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	0.5 mg/ml
Purity	Affinity column purified
Storage	Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRM056 human CD45 [HI30] monoclonal antibody IHC image

Human tonsil tissue



Immunoperoxidase staining of formalin fixed, paraffin-embedded human lymph node and human tonsil tissue.

