

Catalog Number: IRQ6002

Anti-Human CD4 [SP35] monoclonal antibody


Package size: 25, 100 μ l
Store at: -20°C
MW (kDa):

Overview

Product Name	Human CD4 [SP35] monoclonal antibody
Product Number	IRQ6002
Clonality	Monoclonal
Host	Rabbit
Species Reactivity	Human
Recommended Applications Dilutions	Immunofluorescence 1:100 – 1:200 Immunohistochemistry (Paraffin) 1:50 - 1:100
Storage Buffer	Phosphate-buffered solution, pH 7.2, 50% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	0.5 mg/ml
Purity	Affinity column purified
Storage	Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRQ6002 Human CD4 [SP35] monoclonal antibody IHC image

Human tonsil tissue



Immunoperoxidase staining of formalin fixed, paraffin-embedded human tonsil tissue.

